

Growing Tiny Tot)



The GTT (Growing Tiny Tot)
paediatric wheelchair has been
designed to be a child's first
active wheelchair. The unique
design promotes optimal
pushing for the user
allowing them to move
with increased
independence.



USING THE GTT'S WIDE RANGE OF SEAT SIZE OPTIONS AND POSTURE ACCESSORIES, STYLITE DESIGNS WILL SET UP THE CHAIR TO MEET THE USER'S REQUIREMENTS

THE GTT[™] QUICKLY DISASSEMBLES FOR EASE OF TRANSPORT

WITH OVER 20cm OF SEAT WIDTH AND DEPTH ADJUSTMENT. THE GTT™WILL ALLOW FOR ALL YOUR GROWTH REQUIREMENTS

STANDARD PRODUCT FEATURES

- Solid Backrest
- Cushion (5cm thick)
- 2 point pelvic positioning strap
- Choice of four colours (Sunshine Yellow, Fire Engine Red, Fuchsia & Blue)
- · Choice of colourful spoke guards
- · Height adjustable push handles
- · Quick release back wheels
- Flashing front castors
- Detachable backrest
- Angle & height adjustable footrest
- Swing under anti-tippers
- Skirt Guards



CRASH TEST DATA

The Stylite Designs Ltd GTTii Paediatric wheelchair has been crash tested according to the ISO standard 7176/19 Wheelchair Test.

The GTTii was tested with seat configuration (17"x16"). A Hybrid III 5%ile (54kg) crash test dummy was used for the test.

- Wheelchair restraints used: UNWIN
- Occupant restraint used: UNWIN

SPECIFICATIONS	
Weight Limit:	54kg
Seat Width:	20 - 45 cm (2.5 cm increments)
Seat Depth:	15 – 47.5 cm
Back Height:	20 – 45 cm
Leg Length:	15 – 42 cm
Frame Material:	PTG700 High Tensile Steel & Aluminium
Rear Wheels:	18, 20, 22"
Front Wheels:	3,4,5"





PERSONALISED SPOKE GUARDS ALSO AVAILABLE

PLEASE CONTACT US FOR DETAILS.