

PELVI.LOC

2D - Active (2DAS)

For moulded seats

Webbing straps at sides

Code & Size Guide:

PL-2DAS-M (hip width 30-45 cm)

PL-2DAS-L (hip width 36-46 cm)



3DS

For moulded seats

Webbing straps at sides

Code & Size Guide:

PL-3DS-XXXS (hip width 8-16 cm)

PL-3DS-XXS (hip width 12-20 cm)

PL-3DS-XS (hip width 16-24 cm)

PL-3DS-S (hip width 20-32 cm)

PL-3DS-M (hip width 30-40 cm)

PL-3DS-L (hip width 36-46 cm)



3D-Active (3DAS)

For moulded seats

Webbing straps at sides

Code & Size Guide:

PL-3DAS-XS (hip width 16-24 cm)

PL-3DAS-S (hip width 20-32 cm)

PL-3DAS-M hip width 30-40 cm)

PL-3DAS-L (hip width 36-46 cm)



Pelvic retraction seating system

Code & Size Guide:

Pelvi-Contour- seating system Size M
(Hip width 36cm to 42cm)

Pelvi-Contour- seating system Size L
(Hip width 36cm to 46cm)



Neoprene Harness Classic

Code:

PL—W—S

PL—W—M

PL—W—L

PL—W—XL



Ankle Strap

Code & Size Guide:

PL-KG-KS-S - dia. 17-23 cm

PL-KG-KS-M—dia. 20-26 cm

PL-KG-KS-L— dia. 24-30 cm



For more product info please go to: <http://www.a2j-intl.com/en/products-reha/pelvi-loc-system.html>



PELVI.LOC

2D - Active (2DA)

For moulded seats

Plastic straps at sides

Code & Size Guide:

PL-2DA-M (hip width 30-45 cm)

PL-2DA-L (hip width 36-46 cm)



3D – Basic (3DB)

For moulded seats

Plastic straps at sides

Code & Size Guide:

PL-3DB-XS (hip width 16-24 cm)

PL-3DB-S (hip width 20-32 cm)

PL-3DB-M (hip width 30-40 cm)

PL-3DB-L (hip width 36-46 cm)



3D-Active (3DA)

For moulded seats

Plastic straps at sides

Code & Size Guide:

PL-3DA-XS (hip width 16-24 cm)

PL-3DA-S (hip width 20-32 cm)

PL-3DA-M hip width 30-40 cm)

PL-3DA-L (hip width 36-46 cm)



2D Basic (2DB)

Plastic straps at sides

Code & Size Guide:

PL-2DB-XS (hip width 16-24 cm)

PL-2DB-S (hip width 20-32 cm)

PL-2DB-M (hip width 30-40 cm)

PL-2DB-L (hip width 36-46 cm)



Clamping Clips

Required for:

2DA, 3DB, 3DA & 2DB

Code & Size Guide:

RS – 20 (tube dia. 20mm)

RS – 22 (tube dia. 22mm)

RS – 25 (tube dia. 20mm)



PB-BL

Required for:

2DAS, 3DS, & 3DAS



For more product info please go to: <http://www.a2j-intl.com/en/products-reha/pelvi-loc-system.html>

Stylite Designs, Unit 37, Eastage Drive, Little Island, Co. Cork, Ireland

Tel. + 353 (0)21 4355821 Email. info@stylitedesigns.com